

# INNOVATION IN OB/GYN CLINIC

Compared to the alternatives, freezing of condylomata acuminata (genital warts) and precancerous cells on the cervix is much quicker.

Using CryoPen for this treatment is not only very simple but it also brings excellent result in situations when there is repeated bleeding or at the presence of cervical polyps, when HPV is positive, or in the cases where patients suffers from repeated vaginitis..\*

The treatment process goes easy and patient satisfaction is significant.



**CryoPen**®

## SO, WHY HESITATE ?

CryoPen will vitalize your treatment offer with a brilliant, rapid stress free solution that is incredibly accurate and practical.

*"I see this technology as the new standard by which other products for condyloma destruction will be judged." \*\**

These are the features :

- Most convenient ergonomic design from the viewpoint of the doctor
- Clear advantages the cases of multiple viral infections
- No preparation required, no side effects and no follow-up care necessary

### CryoPen | xgyn

With 7 applicators

						
Blue dot ● Ø 1-3mm	White dot ○ Ø 2-6mm	Green dot ● Ø 4-8mm	Yellow dot ● Ø 7-20mm	60mm White ○ Ø 2-6mm	60mm Green ● Ø 4-8mm	120mm Green ● Ø 4-8mm
Power : █	Power : ██	Power : ███	Power : ████	Power : ██	Power : ███	Power : ███

\*Treatment review by Dr. Hae-Lung Sung, Gynecologist, Seoul, Korea.

\*\*Product review by Dr. James Greenberg, chief division gynecology Brigham & Women's Faulkner Hospital and associate Professor, Harvard medical School.